



# INDIAN INSTITUTE OF TECHNOLOGY HYDERABAD

BIOMEDICAL IMAGING

# CT IMAGING

Abhishek Amit Raje

February 2024

## 1 *Sinogram of a Shepp Logan Phantom*

- A sinogram is a graphical representation of the intensity values obtained from a set of projections of an object. It is the projected sum of the 1D Fourier transforms taken about the image.
- The sinogram is represented as a 2D plot where the x-axis corresponds to the projection angle  $\theta$  and the y-axis corresponds to the radial distance  $\rho$ . Each point in the sinogram represents the integral of the the image .

```
1 import numpy as np
2 import matplotlib.pyplot as plt
3 from skimage.data import shepp_logan_phantom
4 from skimage.transform import radon
5
6 def generate_sinogram_with_theta(phantom, theta, cmap='
  viridis'):
7     sinogram = radon(phantom, theta=theta, circle=True)
8
9     plt.subplot(121)
10    plt.imshow(phantom, cmap='gray')
11    plt.title('Shepp-Logan Phantom')
12
13    plt.subplot(122)
14    plt.imshow(sinogram, cmap=cmap, aspect='auto', extent=[
  min(theta), max(theta), 0, sinogram.shape[0]])
15    plt.colorbar(label='Intensity')
16    plt.title('Sinogram' )
17
18    plt.show()
19
20 # Generate Shepp-Logan phantom
21 phantom = shepp_logan_phantom()
22
23 # Specify the sampling angle theta (in degrees) with a
  sampling rate of 0.05 degrees
24 theta_values = np.arange(0, 180,0.05)
25
26 # Generate sinogram with specified theta values
27 generate_sinogram_with_theta(phantom, theta_values)
```

Listing 1: Python code for generating sinogram




Figure 1: Sinogram at a sampling rate  $\theta = 0.5$  degree

## 2 *Filtered Back Projection*

- Filtered back projection is a reconstruction technique to obtain the image after taking the radon transform
- It is used to convert the sinogram back into the 2-D image that the CT scan has taken.
- A filter is applied to the sinogram data to bring attention to other prominent features(frequency) in the image.

```

1 import numpy as np
2 import matplotlib.pyplot as plt
3 from scipy.ndimage import convolve1d
4 from skimage.data import shepp_logan_phantom
5 from skimage.transform import radon
6
7 def generate_sinogram_with_theta(phantom, theta, cmap='
  viridis'):
8     sinogram = radon(phantom, theta=theta, circle=True)
9
10    # Apply Mexican Hat filter
11    mexican_hat_kernel = np.array([-1, 2, -1])
12    filtered_sinogram = convolve1d(sinogram,
13    mexican_hat_kernel, axis=0, mode='constant', cval=0)
14
15    fig, axes = plt.subplots(1, 3, figsize=(15, 5)) # 1 row
16    , 3 columns

```

```

15 axes[0].imshow(phantom, cmap='gray')
16 axes[0].set_title('Shepp-Logan Phantom')
17
18
19 im1 = axes[1].imshow(sinogram, cmap=cmap, aspect='auto',
20 extent=[min(theta), max(theta), 0, sinogram.shape[0]])
21 axes[1].set_title('Original Sinogram')
22 plt.colorbar(im1, ax=axes[1], label='Intensity') # Add
23 colorbar to the first subplot
24
25 im2 = axes[2].imshow(filtered_sinogram, cmap=cmap,
26 aspect='auto', extent=[min(theta), max(theta), 0,
27 filtered_sinogram.shape[0]])
28 axes[2].set_title('Sinogram with Mexican Hat Filter')
29 plt.colorbar(im2, ax=axes[2], label='Intensity') # Add
30 colorbar to the second subplot
31
32 plt.subplots_adjust(wspace=0.5) # Adjust the width
33 space between subplots
34
35 plt.show()
36
37 # Generate Shepp-Logan phantom
38 phantom = shepp_logan_phantom()
39
40 # Specify the sampling angle theta (in degrees) with a
41 sampling rate of 0.05 degrees
42 theta_values = np.arange(0, 180, 0.05)
43
44 # Generate sinogram with specified theta values and apply
45 Mexican Hat filter
46 generate_sinogram_with_theta(phantom, theta_values)

```

Listing 2: Python code for generating sinogram with Mexican Hat filter




Figure 2: Sinogram with Mexican Hat Filter

### Comparing Properties

- As Observed the Mexican Hat filter enhances certain frequency components in the sinogram while suppressing others.
- The sinogram with the Mexican Hat filter is has sharper features
- The filtering process reduces noise in the sinogram, leading to a cleaner representation of the object's characteristics

]

### 3 *Reconstruction of Image*

```
1 import numpy as np
2 import matplotlib.pyplot as plt
3 from skimage.data import shepp_logan_phantom
4 from skimage.transform import radon, iradon
5
6 def generate_sinogram_with_theta(phantom, theta, cmap='
    viridis'):
7     sinogram = radon(phantom, theta=theta, circle=True)
8     plt.subplot(121)
9     plt.imshow(phantom, cmap='gray')
10    plt.title('Shepp-Logan Phantom')
11
12    plt.subplot(122)
13    plt.imshow(sinogram, cmap=cmap, aspect='auto', extent=[
14        min(theta), max(theta), 0, sinogram.shape[0]])
15    plt.colorbar(label='Intensity')
16    plt.title('Sinogram' )
17
18    plt.show()
19
20    return sinogram
21
22 def reconstruct_image_from_sinogram(sinogram, theta):
23     reconstructed_image = iradon(sinogram, theta=theta,
24     circle=True)
25
26     plt.imshow(reconstructed_image, cmap='gray')
27     plt.title('Reconstructed Image')
28     plt.show()
29
30 # Generate Shepp-Logan phantom
31 phantom = shepp_logan_phantom()
32
33 # Specify the sampling angle theta (in degrees) with a
34 # sampling rate of 0.05 degrees
35 theta_values = np.arange(0, 180, 5)
36
37 # Generate sinogram with specified theta values
38 sinogram = generate_sinogram_with_theta(phantom,
39     theta_values)
40
41 # Reconstruct the image from the sinogram
42 reconstruct_image_from_sinogram(sinogram, theta_values)
```

Listing 3: Reconstructing Sinogram using Inverse Radon Transform




Figure 3: Reconstruction with Filter




Figure 4: Reconstruction at sampling angle  $\theta = 0.5$  degree i.e 360 projections

## 4 *Sinogram using 360-degree projection*

### Effect of Sampling Rate on Image Reconstruction

- Insufficient sampling can lead to decreased image reconstruction quality
- When the angle of sampling is too sparse, multiple white streaks appear in the image due to lack of data in the unknown region
- Increasing the number of projections (higher sampling rate) will improve the resolution and characteristics of the reconstructed image.

## 5 *Reconstruction Using Central Slice Theorem*



```

1 import numpy as np
2 import matplotlib.pyplot as plt
3 from skimage.data import shepp_logan_phantom
4 from skimage.transform import radon, iradon
5
6 # Generate Shepp-Logan phantom
7 phantom = shepp_logan_phantom()
8
9 # Perform Radon transform to obtain sinogram
10 sinogram = radon(phantom)
11
12 # Apply Central Slice Theorem by taking the inverse Radon
    transform
13 reconstructed_image = iradon(sinogram)
14
15 # Display the results
16 plt.figure(figsize=(10, 4))
17
18 plt.subplot(131)
19 plt.imshow(phantom, cmap='gray')
20 plt.title('Shepp-Logan Phantom')
21
22 plt.subplot(132)
23 plt.imshow(sinogram, cmap='viridis', aspect='auto', extent
    =[0, 180, 0, sinogram.shape[0]])
24 plt.title('Sinogram')
25
26 plt.subplot(133)
27 plt.imshow(reconstructed_image, cmap='gray')
28 plt.title('Reconstructed Image')
29
30 plt.show()

```

Listing 4: Reconstructing Sinogram using Central Slice Theorem

- the Central Slice Theorem states that the Fourier transform of a 1D projection of an object is equivalent to a central line in the 2D Fourier transform of the object itself.




Figure 5: Reconstruction using central Slice Theorem

## 6 *Reconstruction After Filtering*

```

1 import numpy as np
2 import matplotlib.pyplot as plt
3 from skimage import io, color, transform, img_as_float
4 from skimage.metrics import mean_squared_error
5
6 # Read the image
7
8 # Convert the image to grayscale if it is RGB
9 if image.ndim == 3:
10     image = color.rgb2gray(image)
11
12 # Convert the image to double
13 image = img_as_float(image)
14
15 # Display the size of the matrix
16 print('Size of the image:')
17 print(image.shape)
18
19 # Define the filters
20 filters = ["ram-lak", "shepp-logan", "cosine", "hamming"]
21
22 # Initialize an array to store the mean squared errors
23 mse = np.zeros(len(filters))
24
25 # Loop over each filter
26 for i, filter_type in enumerate(filters):
27     # Perform the radon transform
28     sinogram = transform.radon(image, theta=np.arange(180))
29
30     # Apply the filter
31     filtered_sinogram = transform.iradon(sinogram, filter=

```

```

filter_type)
32
33     # Convert the reconstructed image to double
34     filtered_sinogram = img_as_float(filtered_sinogram)
35
36     # Compute the mean squared error
37     mse[i] = mean_squared_error(image, filtered_sinogram)
38
39     # Display the filter name and corresponding mean squared
40     # error
41     print(f"{filters[i].capitalize()} Filter MSE: {mse[i]:.4f}")
42
43     # Display the images
44     plt.figure(figsize=(12, 4))
45     plt.subplot(131)
46     plt.imshow(image, cmap='gray')
47     plt.title('Original Image')
48
49     plt.subplot(132)
50     plt.imshow(filtered_sinogram, cmap='gray')
51     plt.title(f'Reconstructed Image ({filter_type.capitalize()} Filter)\nMSE: {mse[i]:.4f}')
52
53     plt.subplot(133)
54     plt.plot(np.arange(180), sinogram.mean(axis=0), label='Sinogram')
55     plt.legend()
56     plt.title('Sinogram')
57     plt.show()

```

Listing 5: Reconstructing Sinogram After Filtering

Mean Square Error of Filters:

- Filter: Ram-Lak, Mean Squared Error: 0.0006536
- Filter: Shepp-Logan, Mean Squared Error: 0.0006555
- Filter: Cosine, Mean Squared Error: 0.0007248
- Filter: Hamming, Mean Squared Error: 0.0008254

We can analyze the filtering accuracy of the filter based on the the mean square error

## **7    *Limitations of projection radiography***

- Projection radiography fails to provide details about the 3 dimensional structure as it provides a 2 dimensional image based on the central slice theorem.
- projection radiography provides static images and not suitable for real time imaging
- There are multiple possibilities of distortion of the image due to external factors like patient movement, equipment malfunction.

# X-RAY IMAGING

Abhishek Amit Raje

February 2024

## 1 *Factor Affecting SNR*

Signal-to-Noise Ratio (SNR) in X-ray imaging is influenced by various factors

- **Scatter Radiation:** Scatter radiation contributes to noise in the image
- **Detector Efficiency:** The efficiency of the detector in converting X-ray photons
- **Patient Thickness and Density:** The composition and thickness of the patient's body influence the reduction in intensity of X-rays

If the SNR doubles then dose becomes four times

## 2 *Beam Hardening*

- Beam hardening is a phenomenon in X-ray imaging that occurs when a X-ray beam, which consists of a range of energy levels, passes through an object or a patient. As the X-ray beam passes through different materials of varying thickness and composition
- In this process lower-energy X-ray photons are preferentially absorbed, leading to an increase in the average energy of the remaining X-ray beam.

## 3 *Aliasing, bandwidth limiting, Nyquist condition*

- **Aliasing:** Aliasing occurs when a signal is under sampled, causing incomplete reconstruction of the signal waveform.
- **Bandwidth limiting** Bandwidth limiting restricts the range of frequencies in a signal, between the bandpass range for avoiding interference of irrelevant frequencies.
- **Nyquist condition** The Nyquist condition sampling theorem states that the sampling rate must be at least twice the maximum frequency in a signal to prevent aliasing. This condition is  $f_{sampling} > 2f_{max}$

## 4 *Percentage of X-rays are transmitted through the chest*

We know percentage of X-rays transmitted through the chest, the exponential attenuation formula:

$$I = I_0 e^{-ux} \quad (1)$$

where:

- I is the intensity of the transmitted X-rays

- $I_o$  is the initial intensity of the incident X-rays
- $u$  is the linear attenuation coefficient
- $x$  is the thickness of the material.

The linear attenuation coefficient  $u$  can be calculated using the half-value layer (HVL) ( $H$ ) and the formula

$$u = \frac{\ln 2}{H} \quad (2)$$

Calculate  $u$  for muscle and bone

$$\mu_{\text{muscle}} = \frac{0.693}{3.5 \text{ cm}} \approx 0.198 \text{ cm}^{-1}$$

$$\mu_{\text{bone}} = \frac{0.693}{1.8 \text{ cm}} \approx 0.385 \text{ cm}^{-1}$$

$$I_{\text{muscle}} = I_0 \cdot e^{-\mu_{\text{muscle}} \cdot x_{\text{muscle}}} = I_0 \cdot e^{-0.198 \cdot 16}$$

$$I_{\text{bone}} = I_0 \cdot e^{-\mu_{\text{bone}} \cdot x_{\text{bone}}} = I_0 \cdot e^{-0.385 \cdot 4}$$

$$I_{\text{total}} = I_{\text{muscle}} \cdot I_{\text{bone}}$$

Percentage of X-Ray Transmitted

$$e^{-0.198} * e^{-0.385} = e^{-0.583} = 0.558$$

## 5 *Difference between X-Ray Images*

- Higher Effective Energy (140 keV): X-rays with higher energy penetrate tissues more effectively. This results in lower contrast in the image because the X-rays making it challenging to distinguish between different
- As observed the high energy X-Ray is completely
- Higher Effective Energy (140 keV): The X-rays with higher energy can penetrate the body more deeply, reaching the detector with less absorption of the X-Ray
- Lower Effective Energy (50 keV): Lower energy X-rays are more likely to be absorbed by tissues, leading to greater absorption by tissues.

## 6 *Opinion on getting Xray imaging for heart related diseases.*

Pros	Cons
<ul style="list-style-type: none"><li>• Diagnostic Tool: X-ray imaging can provide valuable information about the structure of the heart.</li><li>• Quick and Non-Invasive: X-rays are relatively quick and non-invasive compared to some other imaging modalities.</li><li>• Widely Available: X-ray facilities are widely available, making it a convenient option for initial screening.</li></ul>	<ul style="list-style-type: none"><li>• Radiation Exposure: X-rays involve exposure to ionizing radiation, which carries potential risks.</li><li>• Limited Soft Tissue Detail: X-rays are better at visualizing bones and dense structures but have limitations in capturing detailed images of soft tissues.</li><li>• May Require Contrast Agents: In some cases, contrast agents may be used, and individuals may have allergies or adverse reactions to these agents.</li></ul>



## 7 *Issue with X-Ray Machine*

The diagram incorrectly labeled the Multitap AC transformer as a "Multiple AC transformer" and the Timing circuit as a "Counting circuit."


The operation of an X-ray machine involves various components:

- Multitap AC transformer: This component adjusts for incoming line variations by selecting different taps to compensate. The number of outputs, or taps, can range from 2 to many, allowing for the selection of a higher or lower voltage tap based on the required X-ray exposure intensity.
- X-ray tube filament transformer: This transformer converts the AC line to supply power for heating the cathode filament. By selecting different taps, the filament heat can be adjusted, subsequently altering the X-ray tube current and total energy delivered to the patient.
- X-ray tube high voltage transformer and bridge rectifier: This unit transforms the AC line into high DC voltage, essential for accelerating electrons from the cathode to the anode. The high DC voltage can be chosen using taps.
- Timing circuit: Contrary to the mislabeling as a "Counting circuit," the Timing circuit controls the timing aspects of X-ray exposure. It includes an electronic counter that regulates the turn-on, turn-off, and duration of X-ray exposure to the patient.

## 8 *Fourier Transforms and Spectra*

a.  $\cos^2(k_0 z)$

$$\mathcal{F}\{\cos^2(k_0 z)\} = \frac{1}{2} [\delta(\omega - k_0) + \delta(\omega + k_0)]$$



b.  $\sin^3(k_0 z)$

$$\mathcal{F}\{\sin^3(k_0 z)\} = \frac{3}{4i} [\delta(\omega - 3k_0) - \delta(\omega - k_0)]$$

